

STOPPER FOR SPIRITS

Stopper for spirits with body made of natural or microagglomerated cork, with wooden or plastic head



| CHARACTERISTICS | PARAMETERS | TOLERANCES |
|-----------------------------------|---|-----------------------|
| Dimensionals | Height of the cork body | Nom. value +/- 0,3 mm |
| | Diameter of the cork body | Nom. value +/- 0,3 mm |
| | Height of the stud | Nom. value +/- 0,2 mm |
| | Diameter of the stud | Nom. value +/- 0,2 mm |
| Chemical / Physical | Humidity | Min: 3% Max: 8% |
| | Peroxides | < 0,2 mg/stopper |
| Sensory | Organoleptict test | <u>≤</u> 2% |
| Features | Dust transfer | < 1,5 mg/stopper |
| Components and treatments | Steam treatment for the reduction of TCA and of other volatile compounds | |
| | Washing: with a solution of H ₂ O ₂ Drying: forced warm air for the elimination of peroxides residuals | |
| | Silicone coating fit to food contact | |
| References to rules and standards | Rules of the "Disciplinare sulla produzione ed utilizzo del tappo di sughero ad uso enologico" ed 2011; sampling procedures according to UNI ISO 2859. In compliance with European and international directives and rules as regarding the food contact. | |
| Shelf-life | Best before 1 year from packaging (storage conditions pointed out in the handbook "the cork - technical handbook for the correct use of cork stoppers"). | |